


# PGS-123

## Pulsed Galvanic Stimulator

- ▶ *Single channel (2-lead) analog high volt pulsed galvanic stimulator (HVPG)*
- ▶ *One active two-lead and one dispersive single-lead output receptacles*
- ▶ *Alternating or continuous treatment mode switch*
- ▶ *Positive or negative polarity switch*
- ▶ *Adjustable frequency (2 - 100 Hz)*
- ▶ *AC wall adapter or 9V battery power supply*
- ▶ *LED output indicator*
- ▶ *Durable metal belt clip*
- ▶ *1-year warranty*



*Comes with 1 each 4" x 7" self-adhesive dispersive electrode, 4 self-adhering 1.5" x 2.5" oval active electrodes (electrodes not shown), 1 each 36" single-lead dispersive lead wire, 1 each 36" 2-lead active lead wire, 1 9-volt heavy duty battery, AC wall adapter, carrying case, and instruction manual*

### Technical Specifications

**Channel** Single, with 1 dispersive and 2 active leads  
**Output Voltage** Max 350 V DC  
**Waveform** Twin peak monophasic  
**Pulse Rate** 2 - 100 Hz  
**Pulse Width** 5  $\mu$ S (fixed)  
**Output Mode** Alternating or continuous  
**Polarity** Positive or negative  
**Power Supply** 9-volt battery or AC wall adapter  
**Size** 6.3 x 9.5 x 3 cm  
**Weight** 116 grams

Distributed by

**TENS Rx**  
Your Real Time Healthcare Solution

Phone: 877-372-4441 | Fax: 877-372-5551